



A NEW AND IMPROVED CHEMILUMINESCENT SUBSTRATE

(Mark Sandison¹, Renuka De Silva¹, Gamachu Melkamu², Patrick Kilmartin², Jaya Koti²)
¹Lumigen Inc., Southfield, MI 48033, ²Beckman Coulter, Chaska, MN 55318

BACKGROUND

Beckman Coulter is developing a new immunoassay system that will run current Access immunoassays as well as additional new menu. Goals for this new system include improved turn-around-times for all assays, thereby meeting STAT test requirements while improving overall platform throughput. A key component of the new system is a new chemiluminescent substrate employed to generate the light signal response. This new substrate is composed of a buffered surfactant enhancer system supporting an alkaline phosphatase-sensitive acridan. When the acridan is triggered in-situ, it forms a dioxetanone which immediately decomposes and emits light.

Lumi-Phos 530 (also known as LP-530) has desirable sensitivity, background luminescence and open bottle stability, but needs 6.3 minutes for signal generation on automated immunoassay systems. The new substrate formulation was optimized for immunoassay specificity, compatibility, sensitivity and is suitable for use with all forms of ALP employed by Access assays with a much shorter time to signal generation.

Comparison of the Lumi-Phos 530 substrate to the new chemiluminescent substrate LumiFAST was done on a immunoassay prototype analyzer to understand the performance characteristics.




Figure 1 illustrates the active component and mechanism for chemiluminescence of Lumi-Phos 530 and LumiFAST substrates




Figure 2 illustrates the emission spectrum of both substrates with an emission maxima around 430 nm for LumiFAST (blue) and 530 nm for Lumiphos-530 (red)




Figure 3 illustrates the signal generation for bovine ALP with both substrates. Peak intensity was seen with LumiFAST within 60 seconds after injection. Results obtained on a BMG Plate Reader over 6 minutes

METHODS

LumiFAST formulation was optimized to work with Access immunoassays. Luminometer read time was assessed by determining the change in relative light unit (RLU) signal over 9 to 72 seconds using an ALP-based enzyme test method and several commercialized Access immunoassays. Improved signal-to-noise performance was demonstrated by comparing calibration curves from several immunoassays generated using Lumi-Phos 530 and the new chemiluminescent substrate LumiFAST. The impact of non-specific signal from endogenous ALP was determined by assessing a panel of patient samples previously identified to contain these interferents, using assays tested with both substrates.

“Assay results shown were generated using immunoassay prototype systems and may not represent final product claims”

RESULTS

Luminometer read time is approximately 5 minutes shorter for the new substrate than for Lumi-Phos 530. Three- to six-fold increases in signal-to-noise performance were demonstrated across the Access immunoassays. Samples with known high endogenous ALP activity displayed greater than 50% reduction in spurious elevations (fliers) when using the new chemiluminescent substrate as compared to the values observed with the same samples using Lumi-Phos 530.

Lumi-Phos 530 on Access 2	LumiFAST on Immunoassay Prototype analyzer
6.3 Minutes signal generation	1 Minute signal generation
Minimum 18 hour room temperature equilibration before use	No room temperature equilibration before use

Calibration Curve Signals





Figure 4 and 5 illustrates the signal generated with hsTnI and TSH calibration curves for both LumiFAST (blue) and Lumi-Phos 530 (red)

Improvement in assay sensitivity (signal-to-noise) with LumiFAST compared to Lumi-Phos 530




Figure 6 illustrates the signal-to-noise improvement across the current immunoassay menu. A 2.5 fold median increase in signal-to-noise was seen with LumiFAST compared to Lumi-Phos 530

Improved Sensitivity with LumiFAST on prototype immunoassay analyzer

TSH Detection Capability (mIU/L)

LumiFast	LoB	LoD	LoQ 10% CV	LoQ 20% CV
Mean	0.00035	0.00074	0.0034	0.00130
Std Dev	0.00008	0.00018	0.0009	0.00031

Lumi-Phos 530	LoB	LoD	LoQ 10% CV	LoQ 20% CV
Mean	0.00073	0.00157	0.0079	0.0028
Std Dev	0.00017	0.00038	0.0029	0.0008




Figure 7 illustrates differences and distinction in RLU between low concentration TSH samples and zero calibrator with both substrates

hsTnI Detection Capability (pg/ml)

LumiFAST	LoB	LoD	LoQ 10% CV	LoQ 20% CV
Mean	0.25	0.32	0.43	0.19
Std Dev	0.10	0.11	0.12	0.05

Lumi-Phos 530	LoB	LoD	LoQ 10% CV	LoQ 20% CV
Mean	0.50	0.51	1.91	0.78
Std Dev	0.31	0.23	0.44	0.16







Figure 8A-D illustrates the Limit of Blank and Limit of Detection for hsTnI with both substrates on the immunoassay prototype analyzer. Better detection capability was seen with LumiFAST





Figure 9A and 9B illustrates differences in LoQ (20% CV) for hsTnI (pg/ml) on immunoassay prototype analyzer

Improved Time to first result with LumiFAST on immunoassay prototype analyzer

Time to first result in minutes

Assay Name	Lumi-Phos 530	LumiFAST
Intact PTH	14	8
Total βHCG	17	11
hsTnI	17	11

Enhanced Assay Specificity with LumiFAST on immunoassay prototype analyzer

An in-house model of incorrect primary tube handling, which leads to neutrophil contamination of plasma samples, is referred to here as the “disturbed tube model”. Neutrophil ALP generates non-specific signal in some immunoassays. False positive results due to this endogenous ALP was evaluated by testing samples subjected to the “disturbed tube model” using normal plasma. HBsAg reactivity (S/CO) with both LumiFAST and Lumi-Phos 530 is shown below

Where (Disturbed == “Yes”)

Signal/Cut-off Distribution of Disturbed Primary Tube by Substrate




Figure 10A and 10B illustrates fewer reactive samples with LumiFAST Substrate compared to Lumi-Phos 530

Disturbed Tube Model- HBsAg

Disturbed	LumiFast		Lumi-Phos 530	
	Non-Reactive	Reactive	Non-Reactive	Reactive
No	74	1	68	7
Yes	65	10	29	46

78% Reduction in false reactivities with LumiFAST compared to Lumi-Phos 530

Patient samples screened for high levels of endogenous ALP were used and were tested for HBsAg reactivity using both Lumi-Phos 530 and LumiFAST substrates

Reactive rate of HBsAg with high endogenous ALP Levels

	Non-Reactive		Reactive		Total
	LumiFAST	Lumi-Phos 530	LumiFAST	Lumi-Phos 530	
LumiFAST	441	30	471		
Lumi-Phos 530	364	111	475		

73% reduction in false reactivities were seen with LumiFAST.

Better Discrimination of Non-Reactive Samples from Cut-off with LumiFAST on immunoassay prototype analyzer

Approximately 500 presumed non-reactive samples were tested using the HIV and HBsAg assays with both Lumi-Phos 530 and LumiFAST chemiluminescent substrates





Figure 11A and 11B illustrates differences in signal to cut-off for HBsAg with LumiFAST and Lumi-Phos 530 for presumed normal samples on immunoassay prototype analyzer






Figure 12A and 12B illustrates differences in signal to cut-off for HIV with LumiFAST and Lumi-Phos 530 for presumed normal samples on immunoassay prototype analyzer



CONCLUSION

The chemiluminescent substrate LumiFAST has been optimized to generate signal rapidly, improve signal-to-noise performance, and reduce non-specific background from endogenous alkaline phosphatase (ALP) in comparison to Lumi-Phos 530. This new substrate presents the opportunity to significantly shorten the time to first result while simultaneously improving assay sensitivity.

- Benefits:
- Shortened Time to first result by ~5 minutes
 - Improved assay sensitivity by reducing signal to noise
 - Improved specificity for assay ALP – Reduced magnitude of falsely elevated signals due to endogenous alkaline phosphatase
 - Improved discrimination of non-reactive and reactive results