

## Chemistry Information Sheet

© 2014 Beckman Coulter, Inc. All rights reserved.

## TBIL Total Bilirubin

REF 442745 (300 tests/cartridge)  
REF 476861 (400 tests/cartridge)

For *In Vitro* Diagnostic Use

Rx Only

### ANNUAL REVIEW

Reviewed by	Date	Reviewed by	Date

## PRINCIPLE

### INTENDED USE

TBIL reagent, when used in conjunction with UniCel<sup>®</sup> DxC 600/800 System(s) and SYNCHRON<sup>®</sup> Systems Bilirubin Calibrator, is intended for quantitative determination of total bilirubin concentration in human serum or plasma.

### CLINICAL SIGNIFICANCE

Bilirubin measurements are used in the diagnosis and treatment of liver, hemolytic, hematological, and metabolic disorders, including hepatitis and gall bladder block.

### METHODOLOGY

TBIL reagent is used to measure the total bilirubin concentration by a timed endpoint Diazo method.<sup>1,2,3</sup> In the reaction, the bilirubin reacts with diazo reagent in the presence of caffeine, benzoate, and acetate as accelerators to form azobilirubin.

The SYNCHRON System(s) automatically proportions the appropriate sample and reagent volumes into a cuvette. The ratio used is one part sample to 35 parts reagent. The system monitors the change in absorbance at 520 nanometers. This change in absorbance is directly proportional to the concentration of TBIL in the sample and is used by the System to calculate and express TBIL concentration.

### CHEMICAL REACTION SCHEME



E015257L.EPS

# SPECIMEN

## TYPE OF SPECIMEN

Biological fluid samples should be collected in the same manner routinely used for any laboratory test.<sup>4</sup> Freshly drawn serum or plasma are the preferred specimens. Acceptable anticoagulants are listed in the PROCEDURAL NOTES section of this chemistry information sheet. Whole blood or urine are not recommended for use as a sample.

## SPECIMEN STORAGE AND STABILITY

1. Tubes of blood are to be kept closed at all times and in a vertical position. It is recommended that the serum or plasma be physically separated from contact with cells within two hours from the time of collection.<sup>5</sup>
2. Separated serum or plasma should not remain at room temperature longer than 8 hours. If assays are not completed within 8 hours, serum or plasma should be stored at +2°C to +8°C. If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples should be frozen at -15°C to -20°C. Frozen samples should be thawed only once. Analyte deterioration may occur in samples that are repeatedly frozen and thawed.<sup>5</sup>
3. Bilirubin is photosensitive. Protect samples from light.

Additional specimen storage and stability conditions as designated by this laboratory:

## SAMPLE VOLUME

The optimum volume, when using a 0.5 mL sample cup, is 0.3 mL of sample. For optimum primary sample tube volumes and minimum volumes, refer to the Primary Tube Sample Template for your system.

## CRITERIA FOR UNACCEPTABLE SPECIMENS

Refer to the PROCEDURAL NOTES section of this chemistry information sheet for information on unacceptable specimens.

Criteria for sample rejection as designated by this laboratory:

## PATIENT PREPARATION

Special instructions for patient preparation as designated by this laboratory:

## SPECIMEN HANDLING

Special instructions for specimen handling as designated by this laboratory:

## REAGENTS

### CONTENTS

Each kit contains the following items:

Two Total Bilirubin Reagent Cartridges (2 x 300 tests) or (2 x 400 tests)

### VOLUMES PER TEST

Sample Volume	8 µL
Total Reagent Volume	280 µL
Cartridge Volumes	
A	255 µL
B	25 µL
C	--

### REACTIVE INGREDIENTS

#### REAGENT CONSTITUENTS

Sodium Benzoate	347 mmol/L
Caffeine	173.9 mmol/L
Sulfanilic acid	27 mmol/L
HCl	50 mmol/L
Sodium Nitrite	0.36 mmol/L
Sodium Acetate	609 mmol/L

Also non-reactive chemicals necessary for optimal system performance.

Avoid skin contact with reagent. Use water to wash reagent from skin.

### GHS HAZARD CLASSIFICATION

Total Bilirubin Reagent  
(Compartment A)

WARNING



H319

Causes serious eye irritation.

P280 Wear protective gloves, protective clothing and eye/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

Sodium Acetate 1 - 10%

Sodium Benzoate 1 - 10%

Total Bilirubin Reagent  
(Compartment B)

DANGER



H314 Causes severe skin burns and eye damage.

H317 May cause an allergic skin reaction.

P280 Wear protective gloves, protective clothing and eye/face protection.

P303+P361+P353 IF ON SKIN (or hair): Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

P362+P364 Take off contaminated clothing and wash it before use

Hydrochloric Acid...% 0.1 - 1%

Sulphanilic Acid 0.1 - 1%

Total Bilirubin Reagent  
(Compartment C)

WARNING



H302 Harmful if swallowed.

H401 Toxic to aquatic life

P264 Wash hands thoroughly after handling.

P273 Avoid release to the environment.

P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

Sodium Nitrite 1 - 10%


## EUROPEAN HAZARD CLASSIFICATION

Total Bilirubin Reagent (Compartment B)	C;R35	
	R35	Causes severe burns.
	S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
	S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
Total Bilirubin Reagent (Compartment C)	T;R25	
	R25	Toxic if swallowed.
	S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).


## MATERIALS NEEDED BUT NOT SUPPLIED WITH REAGENT KIT

SYNCHRON® Systems Bilirubin Calibrator  
 Deionized water (low level calibrator)  
 At least two levels of control material  
 Human serum albumin (azide free)

## REAGENT PREPARATION



For P/N 442745 (2 x 300 tests): Quantitatively transfer 100 microliters (0.1 mL) of the contents from the smallest compartment (C) into the center compartment (B).



For P/N 476861 (2 x 400 tests): Quantitatively transfer 200 µL (0.2 mL) of the contents from the smallest compartment (C) into the center compartment (B).

Replace the cartridge caps and **gently** invert the cartridge several times to ensure adequate mixing. Thorough mixing is necessary for successful calibration.

## ACCEPTABLE REAGENT PERFORMANCE

The acceptability of a reagent is determined by successful calibration and by ensuring that quality control results are within your facility's acceptance criteria.

## REAGENT STORAGE AND STABILITY

TBIL reagent when stored unopened at room temperature will obtain the shelf-life indicated on the cartridge label. Once prepared, the reagent cartridge is stable for 30 days at +2°C to +8°C unless the expiration date is exceeded. DO NOT FREEZE.

Reagent storage location:

## CALIBRATION

### CALIBRATOR REQUIRED

SYNCHRON<sup>®</sup> Systems Bilirubin Calibrator  
Deionized water (low level calibrator)

### CALIBRATOR PREPARATION

No preparation is required.

### CALIBRATOR STORAGE AND STABILITY

If unopened, the SYNCHRON<sup>®</sup> Systems Bilirubin Calibrator may be stored at -15°C to -20°C until the expiration date printed on the calibrator bottle. Opened calibrators that are resealed and stored at +2°C to +8°C are stable for 24 hours unless the expiration date is exceeded.

 CAUTION

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HbsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples to be handled as specified in Centers for Disease Control's Biosafety Level 2 guidelines.<sup>6</sup>

Calibrator storage location:

## CALIBRATION INFORMATION

### NOTICE

Since Total Bilirubin is a calibrated chemistry and also requires "quantitative" reagent preparation it is important to follow proper reagent handling, preparation and storage procedures, especially when utilizing the within-lot calibration feature. Before reporting patient results on successive within-lot cartridges, always analyze and review calibration and quality control data.

1. The system must have a valid calibration curve in memory before control or patient samples can be run.
2. Under typical operating conditions the TBIL reagent cartridge must be calibrated every 14 days and also with certain parts replacements or maintenance procedures, as defined in the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual. This assay has within-lot calibration available. Refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual for information on this feature.
3. For detailed calibration instructions, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.
4. The system will automatically perform checks on the calibration and produce data at the end of calibration. In the event of a failed calibration, the data will be printed with error codes and the system will alert the operator of the failure. For information on error codes, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

## TRACEABILITY

For Traceability information refer to the Calibrator instructions for use.

## QUALITY CONTROL

At least two levels of control material should be analyzed daily. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. More frequent use of controls or the use of additional controls is left to the discretion of the user based on good laboratory practices or laboratory accreditation requirements and applicable laws.

The following controls should be prepared and used in accordance with the package inserts. Discrepant quality control results should be evaluated by your facility.

Table 1.0 Quality Control Material

CONTROL NAME	SAMPLE TYPE	STORAGE

## TESTING PROCEDURE(S)

### NOTICE

When using the within-lot calibration feature it is highly recommended that recovery be confirmed on subsequent cartridge(s) from the same lot number by analyzing quality control material prior to analyzing or reporting any patient results.

1. If necessary prepare reagent as defined in the Reagent Preparation section of this chemistry information sheet and load the reagent onto the system.
2. After reagent load is completed, calibration may be required.
3. Program samples and controls for analysis.
4. After loading samples and controls onto the system, follow the protocols for system operations.

For detailed testing procedures, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

## CALCULATIONS

The SYNCHRON System(s) performs all calculations internally to produce the final reported result. The system will calculate the final result for sample dilutions made by the operator when the dilution factor is entered into the system during sample programming.

## REPORTING RESULTS

Equivalency between the SYNCHRON LX and UniCel DxC 600/800 Systems has been established. Chemistry results between these systems are in agreement and data from representative systems may be shown.

### REFERENCE INTERVALS

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals listed below were taken from literature.<sup>7</sup>

Table 2.0 Reference intervals

INTERVALS	SAMPLE TYPE	CONVENTIONAL UNITS	S.I. UNITS
Literature	Serum or Plasma	0.3 – 1.2 mg/dL	5.1 – 20.5 µmol/L

INTERVALS	SAMPLE TYPE	CONVENTIONAL UNITS	S.I. UNITS
Laboratory			

Refer to References (8,9,10) for guidelines on establishing laboratory-specific reference intervals.

Additional reporting information as designated by this laboratory:

## PROCEDURAL NOTES

### ANTICOAGULANT TEST RESULTS

1. If plasma is the sample of choice, the following anticoagulants were found to be compatible with this method based on a study of 20 healthy volunteers:

Table 3.0 Acceptable Anticoagulants<sup>a</sup>

ANTICOAGULANT	LEVEL TESTED FOR IN VITRO INTERFERENCE	AVERAGE PLASMA-SERUM BIAS (mg/dL)
Sodium Heparin	29 Units/mL	NSI <sup>b</sup>
Lithium Heparin	29 Units/mL	NSI
Ammonium Heparin	29 Units/mL	NSI

a Data shown was collected using SYNCHRON CX Systems. Equivalency between SYNCHRON LX Systems has been established by Deming regression analysis to SYNCHRON CX Systems.

b NSI = No Significant Interference (within  $\pm 0.3$  mg/dL or 6%).

2. The following anticoagulants were found to be incompatible with this method:

Table 4.0 Incompatible Anticoagulants<sup>a</sup>

ANTICOAGULANT	LEVEL TESTED FOR IN VITRO INTERFERENCE	PLASMA-SERUM BIAS (mg/dL) <sup>b</sup>
Sodium Citrate	1.7 mg/mL	$\leq -0.8$
Potassium Oxalate/Sodium Fluoride	4.0 / 5.0 mg/mL	$\leq -0.4$

a Data shown was collected using SYNCHRON CX Systems. Equivalency between SYNCHRON LX Systems has been established by Deming regression analysis to SYNCHRON CX Systems.

b Bias is based on worst case instead of average. Plus (+) or minus (-) signs in this column signify positive or negative bias.

## LIMITATIONS

None identified

## INTERFERENCES

1. The following substances were tested for interference with this methodology:

Table 5.0 Interferences

SUBSTANCE	SOURCE	MAXIMUM LEVEL TESTED	OBSERVED EFFECT <sup>a</sup>
Hemoglobin	RBC hemolysate <sup>b</sup>	100 mg/dL	$\leq +0.63$ mg/dL
Lipemia	Intralipid <sup>c</sup>	200 mg/dL	$\leq -0.24$ mg/dL
Azide <sup>d</sup>	NA <sup>e</sup>	5 mg/dL	$\leq +0.24$ mg/dL
Citrate <sup>d</sup>	NA	900 mg/dL	$\leq \pm 0.20$ mg/dL
Oxalate <sup>d</sup>	NA	1000 mg/dL	$\leq \pm 0.20$ mg/dL
Gentisic Acid <sup>d</sup>	NA	5 mg/dL	$\leq +0.24$ mg/dL
Acetoacetate <sup>d</sup>	NA	0.2 mg/mL 1.08 mg/mL	$\leq +0.7$ mg/dL $\leq +3.7$ mg/dL

a Plus (+) or minus (-) signs in this column signify positive or negative interference.

b Hemolysis may cause falsely elevated results.

c Registered trademarks are the property of their respective owners.

d Data shown was collected using SYNCHRON CX Systems. Equivalency between SYNCHRON LX Systems has been established by Deming regression analysis to SYNCHRON CX Systems.

e NA = Not applicable.

2. Lipemic samples >2+ should be ultra-centrifuged and the analysis performed on the infranate.

3. The Naproxen metabolite, O-desmethylnaproxen, has demonstrated a positive interference with the Jendrassik-Grof method for total Bilirubin measurement.<sup>11</sup>
4. Refer to References (12,13,14) for other interferences caused by drugs, disease and preanalytical variables.

## PERFORMANCE CHARACTERISTICS

### ANALYTIC RANGE

The SYNCHRON System(s) method for the determination of this analyte provides the following analytical ranges:

Table 6.0 Analytical Range

SAMPLE TYPE	CONVENTIONAL UNITS	S.I. UNITS
Serum or Plasma	0.1 – 30.0 mg/dL	1.7 – 513.0 µmol/L

Samples with concentrations outside the analytical range will be reported as "<0.1 mg/dL" ("<1.7 µmol/L") or ">30.0 mg/dL" (">513.0 µmol/L").

Samples reported out as greater than the analytical range may be confirmed by diluting with human serum with a known bilirubin value and reanalyzing. The appropriate dilution factor should be applied to the reported result.

### REPORTABLE RANGE (AS DETERMINED ON SITE):

Table 7.0 Reportable Range

SAMPLE TYPE	CONVENTIONAL UNITS	S.I. UNITS

### SENSITIVITY

Sensitivity is defined as the lowest measurable concentration which can be distinguished from zero with 95% confidence. Sensitivity for TBIL determination is 0.1 mg/dL (1.7 µmol/L).

### EQUIVALENCY

Equivalency was assessed by Deming regression analysis of patient samples to accepted clinical methods.

Serum or plasma (in the range of 0.2 to 28.2 mg/dL):

Y (SYNCHRON LX Systems)	= 0.96X + 0.31
N	= 79
MEAN (SYNCHRON LX Systems)	= 6.75
MEAN (SYNCHRON CX7 DELTA)	= 6.69
CORRELATION COEFFICIENT (r)	= 0.9997

Refer to References (15) for guidelines on performing equivalency testing.

### PRECISION

A properly operating SYNCHRON System(s) should exhibit precision values less than or equal to the following:

Table 8.0 Precision Values

TYPE OF PRECISION	SAMPLE TYPE	1 SD		CHANGEOVER VALUE <sup>a</sup>		% CV
		mg/dL	µmol/L	mg/dL	µmol/L	
Within-run	Serum/Plasma	0.15	2.6	5.0	86.7	3.0
<b>Total</b>	Serum/Plasma	0.22	3.8	5.0	86.7	4.5

a When the mean of the test precision data is less than or equal to the changeover value, compare the test SD to the SD guideline given above to determine the acceptability of the precision testing. When the mean of the test precision data is greater than the changeover value, compare the test % CV to the guideline given above to determine acceptability. Changeover value = (SD guideline/CV guideline) x 100.

Comparative performance data for a SYNCHRON LX<sup>®</sup> System evaluated using the NCCLS Proposed Guideline EP5-T2 appears in the table below.<sup>16</sup> Each laboratory should characterize their own instrument performance for comparison purposes.

Table 9.0 NCCLS EP5-T2 Precision Estimate Method

TYPE OF IMPRECISION	SAMPLE TYPE		No. Systems	No. Data Points <sup>a</sup>	Test Mean Value (mg/dL)	EP5-T2 Calculated Point Estimates	
						SD	%CV
Within-run	Serum	Control 1	1	80	1.7	0.1	6.1
	Serum	Control 2	1	80	5.9	0.1	1.7
	Serum	Control 3	1	80	8.9	0.1	1.2
	Serum	Control 4	1	80	17.5	0.2	1.1
<b>Total</b>	Serum	Control 1	1	80	1.7	0.1	6.1
	Serum	Control 2	1	80	5.9	0.1	1.9
	Serum	Control 3	1	80	8.9	0.1	1.3
	Serum	Control 4	1	80	17.5	0.2	1.23

a The point estimate is based on the pooled data from one system, run for twenty days, two runs per day, two observations per run on an instrument operated and maintained according to the manufacturer's instructions.

**NOTICE**

These degrees of precision and equivalency were obtained in typical testing procedures on a SYNCHRON LX<sup>®</sup> System and are not intended to represent the performance specifications for this reagent.

## ADDITIONAL INFORMATION

For more detailed information on UniCel DxC Systems, refer to the appropriate system manual.

Beckman Coulter, the Beckman Coulter Logo, Synchron, UniCel and DxC are trademarks of Beckman Coulter, Inc and are registered in the USPTO.

### SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

### REVISION HISTORY

#### Revision AF

Revised Interferences section.

**Revision AG**

Updated corporate address; updated OSHA precaution, removed insert reference from content description, and removed EDTA as an Acceptable Anticoagulant claim.

**Revision AH**

Added Reagent Preparation visual aid to the Reagent Preparation section.

**Revision AJ**

Added Revision History.

**Revision AK**

Revised Interferences section.

**Revision AL**

Added new language requirement: Czech, and Korean.

**Revision AM**

Removed references to CX and LX systems as they are discontinued effective 12/2013.

Added Beckman Coulter trademark statement and disclaimer.

**Revision AN**

Added GHS Classification information

## REFERENCES

1. Malloy, H. T., Evelyn, K. A., *J. Biol. Chem.*, 119 481 (1937).
2. Jendrassik, L., Grof, P., *Biochem. Z.*, 297:81 (1937).
3. U.S. Patent No. 4,672,041.
4. Tietz, N. W., "Specimen Collection and Processing; Sources of Biological Variation", *Textbook of Clinical Chemistry*, 5th Edition, W. B. Saunders, Philadelphia, PA (2005).
5. National Committee for Clinical Laboratory Standards, *Procedures for the Handling and Processing of Blood Specimens* Approved Guideline, NCCLS publication H18-A, Villanova, PA (1990).
6. CDC-NIH, *Biosafety in Microbiological and Biomedical Laboratories*, 5th Edition, (Washington, D.C.: U.S. Government Printing Office, 2009). (CDC 21-1112)
7. Tietz, N. W., *Clinical Guide to Laboratory Tests*, 3rd Edition, W. B. Saunders Company, Philadelphia, PA (1995).
8. National Committee for Clinical Laboratory Standards, *How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory* Approved Guideline, NCCLS publication C28-A, Villanova, PA (1995).
9. Tietz, N. W., ed., *Fundamentals of Clinical Chemistry*, 6th Edition, W. B. Saunders, Philadelphia, PA (2007).
10. Henry, J. B., *Clinical Diagnosis and Management by Laboratory Methods*, 22nd Edition, W. B. Saunders Company, Philadelphia, PA (2006).
11. Nafila, A.R., et al, *Spurious hyperbilirubinemia caused by naproxen*, *Clinical Biochemistry* 42 129-131 (2009)
12. Young, D. S., *Effects of Drugs on Clinical Laboratory Tests*, 5th Edition, AACC Press, Washington, D. C. (2000).
13. Friedman, R. B., Young, D. S., *Effects of Disease on Clinical Laboratory Tests*, 4th Edition, AACC Press, Washington, D.C. (2001).
14. Young, D. S., *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 3rd Edition, AACC Press, Washington, D. C. (2007).
15. National Committee for Clinical Laboratory Standards, *Method Comparison and Bias Estimation Using Patient Samples* Approved Guideline, NCCLS publication EP9-A, Villanova, PA (1995).
16. National Committee for Clinical Laboratory Standards, *Precision Performance of Clinical Chemistry Devices* Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992).

EC	REP
----	-----

 Beckman Coulter Ireland, Inc., Mervue Business Park, Mervue, Galway, Ireland 353 91 774068

 Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.